

Weekly Updates

February 19, 2015

Test	Description of Changes	Affected Individuals
Dimension	Please review all of the flip chart Job Aids as they have been updated!!	Bench Techs
Dimension Cal	When making up calibrators, use the volumetric pipettes. They are more precise and will reduce the 'shifting' of QC ranges.	All Staff
Internship/ Training	FYI: We will be taking a break from students until this fall. 😊	All Staff
T Spot TB Test	The labels have changed... There is a tube A with specimen type of 'other' – ignore this label . Then tubes B and C are for dark green tubes, label the specimens with these labels. When putting them on the transfer list, <u>put tubes A, B & C on the list</u> .	All Staff

Procedure Additions, Updates and Changes: