

Protocol	Appendix C: Molecular Work Areas by Room Number					
	Reagent Prep / Mastermix	Extraction / Template Prep	Template Addition	Amplification	Post- Amplification	Additional Handling
	Pre-amplification / Assay Set Up	Extraction / Template Addition Area	Amplification Area	Post-Amplification		
Norovirus Duplex RT-PCR	142	143	140	Micro		
Shiga-toxin <i>E. coli</i> (STEC) RT-PCR	142	Micro	140	Micro		
BAX RT-PCR	N/A	Micro	Micro	Micro		
Salmonella Serotyping BioPlex assay	142	Micro	143	Micro	140	Bioplex 200
BT RT-PCR	116	105A	140	Micro		
TB RT-PCR	142	105B/106/143	140	Micro		
Influenza Virus RT-PCR	142	143	140	Micro		
West Nile Virus RT-PCR	142	107	140	Micro		
Cystic Fibrosis	142	NBS	NBS	140	139	read sealed cards - NBS
SCID	142	NBS	NBS	NBS		Stds prep - 140

Appendix C: Molecular Work Areas by Room Number  
revised February 10, 2015