

HSV-Typing PCR Master Mix QC Form

Run Title: _____

MM lot #: _____

Test: HSV Typing PCR Qiagen_____ or ABI_____

Thermocycler: _____

HSV-1:

(Old Typing Hi+ Ct) - (New HSV Typing Hi+ Ct) = ± 1 Acceptable: ___Yes ___NO

(Old Typing Lo+ Ct) - (New HSV Typing Lo+ Ct) = ± 1 Acceptable: ___Yes ___NO

HSV-2:

(Old Typing Hi+ Ct) - (New HSV Typing Hi+ Ct) = ± 1 Acceptable: ___Yes ___NO

(Old Typing Lo+ Ct) - (New HSV Typing Lo+ Ct) = ± 1 Acceptable: ___Yes ___NO

Actions: If the master mix lot is not acceptable, describe the actions taken below.

Date: _____

Date: _____

Tech: _____

Reviewer: _____