UW Medicine UNIVERSITY OF WASHINGTON MEDICAL CENTER

Irradiator Maintenance Log

YEAR	January	February	March	April	Мау	June	July	August	September	October	November	December
Date (M/DD/YY)												
Inspect door interlock mechanism (Y or N)										<u> </u>		
Inspect power cord (Y or N)												
Clean Canisters (Y or N)												
Clean Canister Holder (Y or N)												
Dust top vents (Y or N)												
Export and Save Irradiator Logs												
Dosimetry Testing (Y or N)		N	N	N	N	N		N	N	N	N	N
Tech ID												

KEY: Y = Yes N = No

Comments:

SUPERVISORY REVIEW:	DATF [,]