

Cell #	Rh-ir	Donor Number	Rh-ir										KELL					DUFFY			KIDD		LEWIS			MNS		P		LUTHERAN	Special Antigen Typing	Test Results
			D	C	E	c	e	f	Cw	V	K	K	Kpa	Kpb	Js ^a	Js ^b	Fy ^a	Fy ^b	Jka	Jkb	Xga	Lea	Leb	S	s	M	N	P ₁	Lu ^a			
1	R1R1	317661	+	+	0	0	+	0	0	0	+	+	/	+	+	0	+	0	0	+	0	0	+	+	0	+	0	+				
2	R2R2	319929	+	0	+	+	0	0	0	0	+	+	/	+	+	+	+	+	+	+	0	+	0	+	+	0	+					
3	rr	319826	0	0	0	0	+	+	+	+	0	0	+	+	+	+	0	+	0	0	+	+	0	+	+	+	0	+				
	Patient Cells																															

Shaded columns indicate those antigens which are destroyed or depressed by enzyme treatment.

* represents "Not Tested" for new donors

Ortho Clinical Diagnostics

Reagent Red Blood Cells
0.8% Surgiscreen®
 ©Ortho-Clinical Diagnostics, Inc. 2010

LOT NO. **VSS960**
 EXP. DATE **2018-01-02**
 CCYY-MM-DD
ANTIGRAM®
 Antigen Profile
 635200621

IMA TEST
000-00-XXXX
12-15-17 MH

Cells: 1 2 3
 1 4+
 2 3+
 3 2+

