

WHAT IS NEW IN CERNER UPGRADE 2018.04.02


Blood Bank Applications

1. Application: PRODUCT HISTORY REVIEW


Already in PROD: You can now unlock a locked product in Product History Review.

Demonstration:


Loaded product in Correct Inventory; this will be now be locked:




Load the locked product in Product History Review




Select Unlock Product from the Task menu.



Enter a reason to unlock the product and a comment, if needed




After OK is clicked the product is no longer locked and the lock icon is no longer displayed



2. Application: RECEIVE PRODUCTS


- a. Currently, in PROD, when ALT + E are pressed, the Edit menu does not open when clicked in any field.

Demonstration (ALT+E were pressed):



- b. After the **Upgrade to 2018.04.02**, the Edit option on the toolbar opens the menu correctly when ALT + E are pressed in Receive Products.

Demonstration: (ALT+E were pressed):



3. Application: RECEIVE PRODUCTS

Occasionally, the bar code value was displayed in the Product Type box
 Example: E2555 FP displays on the product type field as E2555V00.

This has now been corrected. The product type will now display in the Product Type box in Receive Products consistently.

