

BMS/Indiana Univ - Credo Cube packaging reminder

Apheresis Operations <apheresis@celltherapy360.com>

Thu 7/7/2022 9:55 AM

To: Sexton, Steven M <ssexton1@IUHealth.org>; Schwering, Dave T <dschwering@IUHealth.org>

**** EXTERNAL Message From apheresis@celltherapy360.com. DO NOT open attachments or click links from unknown senders or unexpected emails. ****



Dear Steven and Dave,

Please share the following general reminder with all staff responsible for packaging BMS products in the Credo Cube shipper.

The BMS Credo Cube shipper contains a corrugated box in the center compartment of the shipper with packaging/shipping supplies. This box should be removed and discarded before placing the apheresis product directly into the center compartment of the cooling engine. An example of the supply box required to be removed and discarded is shown below.

Please refer to Section 9.5, step 9.5.6. of the BMS MNC Collection Procedure which states to remove all items from the shipping container, including the inner cardboard corrugate.



Please contact your Apheresis Operations account manager at apheresis@celltherapy360.com with questions.

Thank you,

Apheresis Operations

[Apheresis@CellTherapy360.com](mailto:apheresis@CellTherapy360.com)

888-805-4555 Option 2, option 1 (US)

855-999-0170 Option 2, option 1 (Canada)



Apheresis Portal: ct360.com

This message (including any attachments) may contain confidential, proprietary, privileged and/or private information. The information is intended to be for the use of the individual or entity designated above. If you are not the intended recipient of this message, please notify the sender immediately, and delete the message and any attachments. Any disclosure, reproduction, distribution or other use of this message or any attachments by an individual or entity other than the intended recipient is prohibited.