




Template: TMP-00582 v1.0

To: Kite Konnect® Apheresis Collection Module (ACM) Users
From: Kite US Site Qualification (Quality) Teams
Date: 26JAN2024
Subject: Kite Konnect® Apheresis Collection Module (ACM) January 2024 Enhancements

To Whom It May Concern:

Effective 29JAN2024, Kite will upgrade its Kite Konnect® Apheresis Collection Module (ACM) with the following enhancements:

Enhancement	Rationale
1) On the final Apheresis Report, users will now see completed activities from the Verify Shipper Kit Stage , Confirm Handoff Stage , and Packaging Summary Stage .	As these are steps which are completed by staff, having them also listed on the final Apheresis Report creates a more holistic record and summary of the Kite collection.
2) On the final Apheresis Report DIN is currently listed only as the alphanumeric ID. With this update, users will now also see a barcode generated in the DIN field (if a DIN is entered). <div data-bbox="253 1037 945 1108" style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p>DIN to be printed on the Kite Final Product Label</p>  </div>	For sites who reference the listed DIN on the final Apheresis Report for further record keeping; adding a scannable element adds efficiency and decreases the potential for transcription errors.
3) As there are some US facilities with multiple collection locations, two new columns were added to the ACM homepage so that “Institution Name” and “Collection Facility” for each patient is easily seen from the homepage. This will be visible on all three patient tabs. Additionally, all collected patients for an institution, regardless of collection facility, will be visible in the “completed patients” tab. <div data-bbox="253 1404 961 1554" style="border: 1px solid black; padding: 5px; margin-top: 10px;"> <p>Active Patients Completed Patients Cancelled Patients</p> <hr/> <p>Patient Name <input type="text"/> Institution Name <input type="text"/> Collection Facility <input type="text"/> Kite Patient ID Date of Birth</p> </div>	This is being for ease of patient management from the ACM homepage.

This memorandum is for your information only. No further action or training is required as the user instructions remain the same. If you have any questions, please contact your Kite Quality representative.

This is an electronically controlled document. User is responsible for verifying and using the current approved version.

APPROVED