

# COI/C 5.1 Release Scope

## Apheresis and Cell Lab Users

May 4, 2025 Go Live

**Johnson&Johnson**




# Landing Page Patient Selection Table Count Update

## DESCRIPTION

Updating the number of cases that are shown on the landing page for improved functionality and visibility of the Order Details Pop-Up

## IMPACT

All portal users will only see 15 cases per page in the patient selection table.

-  Global
-  All Portal Users
-  Portal



	Order ID	COI ID	Patient Name	Lot Number	Est. Delivery ... 1	Order Journey St..	Status	Annotate
<input type="radio"/>	EU-MSL-10129-BZFM-0809	CLC.TW9111	Shepherd Darryl	QW476T	Mar 05, 2025	Manufacturing	In Process	
<input type="radio"/>	US-MSKCCS-0E51-2363	CLC.21A111	Paulette Jennifer	QW476T	Mar 06, 2025	Manufacturing	Scheduled	
<input type="radio"/>	US-MSKCCS-IWT5-3567	CLC.11A111	Jacqueline Jeremiah	QW476T	Mar 06, 2025	Manufacturing	Scheduled	
<input type="radio"/>	US-MSKCCS-UK4C-0829	CLC.WZ9111	Eileen Glenda	QW476T	Mar 06, 2025	Manufacturing	Scheduled	
<input type="radio"/>	US-MSKCCS-T7MT-5538	CLC.7Z9111	Bryan Munoz	QW476T	Mar 06, 2025	Manufacturing	Scheduled	
<input type="radio"/>	EU-MSL-10129-U7MJ-8195	CLC.GY9111	Rios Guy	852941	Mar 11, 2025	Manufacturing	In Process	
<input type="radio"/>	EU-MSL-10129-LNQ0-3571	CLC.EY9111	Deborah Robyn	QW476T	Mar 11, 2025	Manufacturing	In Process	
<input type="radio"/>	EU-TRT-CH-20233-5YXD-6603	CLC.61A111	Eugene Natalie	QW476T	Mar 12, 2025	Manufacturing	Scheduled	
<input type="radio"/>	EU-TRT-CH-20233-TDVI-6633	CLC.51A111	Nathaniel Denise	QW476T	Mar 12, 2025	Manufacturing	Scheduled	
<input type="radio"/>	EU-MSL-10129-L4TH-5164	CLC.WY9111	Sanders Erin	QW476T	Mar 12, 2025	Manufacturing	In Process	
<input type="radio"/>	EU-MSL-10129-3T81-1062	CLC.UY9111	Autumn Mr.	QW476T	Mar 12, 2025	Manufacturing	In Process	
<input type="radio"/>	EU-MSL-10129-DOZI-2399	CLC.RY9111	Nathaniel Angel	QW476T	Mar 12, 2025	Manufacturing	In Process	
<input type="radio"/>	US-MSKCCS-IP475-47372	CLC.KZ9111	Geralt Rivia	Q20250226	Mar 25, 2025	Manufacturing	In Process	Remanufacturing
598 Results					<div><div>&lt;</div><div>1</div><div>...</div><div>40</div><div>&gt;</div></div>			

# Label Re-Generation


## DESCRIPTION


Operators no longer have to select the "Modify" button and re-submit data that wasn't being modified in order to re-generate and re-print labels.

## IMPACT

Only after generating a label at least once and after completing the double signature an operator can utilize a new "Reprint Label" button at the top of the workflow. Upon clicking the button, an operator will be asked to enter a reason for re-printing and click “Submit”.

 Global

 Apheresis, Local Cell Lab, Central Cell Lab Operators

 Collection Shipment, Cryopreservation, Cryopreservation Shipment

J&J

janssen

English

Shipment to Cryopreservation > Data Collection > Summary

Transfer Product

Shipment Checklist

Review

Reprint Label

Modify

Patient Information: Percy Weasley

Completed

First Name	Middle Name	Last Name	Suffix	Date of Birth	Collection Date
Percy		Weasley		12.05.1986	24.02.2025
Treatment Site	COI ID	Order ID	Therapy/Product		
Universitätsklinikum Leipzig AG, GMP-Labor	OLC-AZ7111	EU-TRTTR232124-00017	Cita-Cel		

You have completed this flow successfully and can now download the summary

# Thank You!