

## Whole Blood QC

All 4 tubes used for QC.

Profiles in table and SOP.



#	Profile	Assigned QC Vials	Expected Results
1	ABORH	Vial 3	A Pos
2	ABO3CELL	Vial 1 Vial 4	AB Pos, IAT Neg B Neg, ABS Pos (anti-D)
3	DONOR	Vial 3 Vial 4	A Pos B Neg
4	ABORH2	Vial 2	O Pos
5	3CELL	Vial 2 Vial 1	ABS Pos (anti-c) ABS Neg
6	Weak D	Surgiscreen 3	D Neg
7	Weak D2	Surgiscreen 1	D Pos
8	QC Buffer Card	Vial 3	Group A
9	Diluent 2/XM QC	Ortho QC C Ortho QC I	Compatible MTS XM Incompatible MTS XM

04/01/2026 RO