

# F° Freezer Temperature Log

Record CURRENT, MIN, and MAX temperatures twice a day. Complete steps 1-4.

**Step 1** Write your initials and the time of day.

Staff Initials	EH																		
Day of Month	example	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26
Time	8:00 am	1:10 pm	8:22 am	3:30 pm															

**Step 2** Write CURRENT, MIN, and MAX temperatures. Circle any temperatures that are in the DANGER Zone. Then go to step 3 (even if CURRENT temp is OK). If ALL temperatures are OK, go to step 4.

CURRENT	example	2.1	0.23	1.60															
MIN		-1.05	-3.06	-2.56															
MAX		6.2	10.35	0.13															

-20s° & lower -10s° -9° -8° -7° -6° -5° -4° -3° -2° -1° 0° 1° 2° 3° 4° 5°  
 These temperatures are OK. Go to Step 4.

6° 7° 8° 9° 10° 11° 12° 13° 14° 15° 16° 17° 18° & higher  
 DANGER Zone  
 Too Warm! Go to Step 3!

**Step 3**

- If any CURRENT, MIN, or MAX temperature is in the Danger Zone (above 5°F):
- Alert your supervisor immediately.
- Press the MEMORY CLEAR/RESET button(s).
- Check the temperatures again in 1 hour. If temps are still in the DANGER Zone, call the VFC Call Center (1-877-243-8832) to report the incident.


Document the date and actions you take:

5/14/14 Freezer temp rechecked, at 9:05am Max temp in Danger Zone. Change Nurse (Denise) notified.  
 SAT/14 freezer temp rechecked @ 9:15a max temp in danger zone ON (notified) needed again @ 10:15a temp 0.13° within normal range.

**Step 4**

- Press the MEMORY CLEAR/RESET button(s) on the thermometer every time you finish logging temperatures. If you press the MODE/ALARM button by mistake, LO and HI will display. Press it again to see MIN and MAX.

At the end of the month file this log and keep it for 3 years.



Example of thermometer display

VACCINES FOR CHILDREN CALIFORNIA



www.eziz.org

Month/Year May / 2014  
 (Days 9-16)  
 Freezer location \_\_\_\_\_  
 VFC PIN \_\_\_\_\_