

Phenocell B 3% Antigen Composition Sheet

Batch No: 1664 197 Expiry Date: 23/12/2019

19-032-2155

Cell No	Reference No	RH							KEL			FY		JK		MNS			P1PK	LE		LU	Additional Typings	Cell No	Results		
		Rh Phen.	D	C	E	c	e	C ^w	K	k	Kp ^a	Fy ^a	Fy ^b	JK ^a	JK ^b	M	N	S		s	P1						Le ^a
1	4135608	R ₁ R ₁	+	+	0	0	+	+	0	+	0	0	+	+	0	+	+	0	+	0	+	0	0	1	2	2	
2	2146231	R ₁ R ₁	+	+	0	0	+	0	+	0	+	0	0	+	+	0	+	0	+	0	+	0	0	2	2	2	
3	2227853	R ₁ R ₁	+	+	0	0	+	0	0	+	0	+	0	+	0	0	+	+	+	+	+	0	0	3	2	2	
4	4398907	R ₂ R ₂	+	0	+	+	0	0	0	+	0	0	+	0	+	0	+	0	+	+	+	0	+	4	3	2	
5	2278162	R ₂ R ₂	+	0	+	+	0	0	0	+	+	0	+	0	+	0	+	0	+	0	+	0	0	5	3	2	
6	2172354	r ^r r	0	+	0	+	+	0	0	+	0	+	0	+	+	+	0	+	+	+	+	0	0	6	2	2	
7	2138870	r ^r r	0	0	+	+	+	0	0	+	0	0	+	0	+	+	+	0	+	0	0	+	0	7	3	2	
8	2269224	rr	0	0	0	+	+	0	+	+	0	0	+	0	+	0	+	+	+	+	0	0	0	8	2	2	
9	2055924	rr	0	0	0	+	+	0	0	+	0	0	+	+	+	+	0	+	0	+	0	0	0	9	2	2	
10	2204055	rr	0	0	0	+	+	0	0	+	0	+	0	+	+	+	0	0	+	+	+	0	0	10	2	2	
11	2117666	rr	0	0	0	+	+	0	0	+	0	+	+	0	+	0	+	0	+	0	0	+	0	11	2	2	
Auto																							Auto	3			
															</												

*Cs^a: A high incidence antigen (Ref: Reid M.E., Lomas-Francis C., Olsson M.L., "The Blood Group Antigen FactsBook", 3rd Edition Pgs: 648-649).

Note: Cells 6, 7 and 8 form the RhD Negative Screening Cell Subset, formulated for antenatal RhD antibody screening in the presence of Prophylactic Anti-D.

This product is manufactured from blood donated by voluntary donors to the Australian Red Cross Blood Service (ARCBS).