

TESTING FOR COLD AGGLUTININS

Appendix 2:

Cold Panel Result Sheet

Patient Name : _____ Patient Ur : _____

Patient Accession : _____ Date : _____ Scientist Initial : _____

Reagent Cell	Lot	Expiry	37°C Sal	RT Sal	4°C Sal	Comment
Cord	Immucelb 665 51203	20/7/2020	0	+/-	0.5	I antigen
Auto			0	0	0.5	
A1	ortho A4522 20/7/20		0	0.5	0.5	Include if patient is blood group A or AB Stronger reaction in A1 cells than A2 May indicate an Anti-A1
A2	DONOR 3937800	28/7/2020	0.5	3	3	
B	ortho A4522 20/7/20		0	2	2	Include if patient is blood group B or AB (I and B antigen)
Screening cell 1	ortho 355812	14/7/2020	0.5	3	4	Group O cells (H and I antigen)
Screening Cell 2	ortho 355812	14/7/2020	0.5	3	4	
Screening Cell 3	ortho 355812	14/7/2020	0.5	3	4	

Of Webb 9/7/19