

PANOCELL-20

Master List

IMMUCOR, INC. Norcross, GA 30071 USA

U.S. License No: 886

Lot No: 49841

Exp. Date: 2017/02/10

NAME _____ NO. _____

INSTITUTION _____

BLOOD GROUP _____ ANTIBODY IDENTITY _____

TECH _____ DATE _____

																							PATIENT'S TEST RESULTS												
SYSTEM	Rh - Hr							Kell					Duffy		Kidd		Lewis		P	MN				Lutheran		Xg	Special Antigen Type								
Donor	D	C	c	E	e	V*	C ^w	K	k	Kp ^a	Kp ^b	Js ^{a*}	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b							Xg ^a		
1	R1wR1 B2599	+	+	0	0	+	0	+	+	0	+	0	+	0	+	+	0	0	0	0	0	+	0	+	0	+	+	+							1
2	R1R1 B6447	+	+	0	0	+	0	0	0	+	0	+	0	+	0	+	0	+	0	+	+	+	0	+	0	+	+	+	Lu:14						2
3	R1R1 B9499	+	+	0	0	+	0	0	0	+	0	+	0	+	+	+	0	+	0	0	0	+	0	+	+	0	+	0							3
4	R1R1 B9498	+	+	0	0	+	0	0	0	+	0	+	0	+	0	+	0	+	+	+	+	0	+	0	+	0	+	0							4
5	R1R1 B9503	+	+	0	0	+	0	0	0	+	0	+	0	+	+	0	+	+	0	0	+	+	0	+	0	0	+	+							5
6	RzR2 A4290	+	W	+	+	0	0	0	0	+	0	+	0	+	0	+	0	0	0	+	+	0	+	0	+	0	+	+							6
7	R2R2 C4651	+	0	+	+	0	0	0	0	+	0	+	0	+	+	+	0	+	0	+	+	+	+	0	+	0	+	0							7
8	R2R2 C5901	+	0	+	+	0	0	0	0	+	0	+	0	+	+	+	+	0	+	+	+	+	+	0	+	+	+	Di(a+)							8
9	R2R2 C5663	+	0	+	+	0	0	0	0	+	0	+	0	+	+	0	0	+	+	0	+	+	0	+	+	+	Bg(a+)							9	
10	R1r R3038	+	+	+	0	+	0	0	0	+	0	+	0	+	0	0	+	0	0	0	+	+	0	+	0	+	+							10	
	D	C	c	E	e	V*	C ^w	K	k	Kp ^a	Kp ^b	Js ^{a*}	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a								
11	r'r E458	0	+	+	0	+	0	0	+	+	0	+	0	+	0	+	0	0	+	0	0	+	0	+	0	+	+							11	
12	r'r F730	0	0	+	+	+	0	0	0	+	0	+	0	+	+	+	+	0	+	0	+	+	+	+	0	+	+								12
13	rr H1322	0	0	+	0	+	0	0	0	+	0	+	0	+	+	0	0	+	+	0	+	+	0	0	+	0	+								13
14	rr G1737	0	0	+	0	+	0	0	+	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+	+								14
15	rr N4310	0	0	+	0	+	0	0	0	+	0	+	0	+	0	+	0	0	+	0	+	0	0	+	0	+	+								15
16	rr N4243	0	0	+	0	+	0	0	0	+	0	+	0	+	+	+	+	0	+	0	+	0	+	0	+	0	+								16
17	rr G1541	0	0	+	0	+	0	0	+	+	0	+	0	+	+	+	+	+	0	+	+	+	0	+	+	+	+								17
18	rr H1049	0	0	+	0	+	0	0	0	+	+	+	0	+	+	0	0	+	0	+	+	+	0	+	+	+	+								18
19	rr N4312	0	0	+	0	+	0	0	0	+	0	+	0	+	+	+	+	+	0	+	+	+	0	+	0	+	+								19
20	Ror D1727	+	0	+	0	+	0	0	0	+	0	+	0	+	+	0	+	+	0	+	+	+	+	+	0	+	+								20
PATIENT'S CELLS																							PATIENT'S CELLS												

In those instances where a patient's serum is known to contain anti-D, it may be desirable to perform antibody screening tests with D- red cells. The panel cells 11,12 and 13 can be used together to form a D- negative antibody screening reagent.

* Indicates those antigens whose presence or absence may have been determined using only a single example of a specific antibody. An antigen designated with a 'w' represents a weakened expression of the antigen that may or may not react with all examples of the corresponding antibody.

REVERSE GROUPING CELLS	A ₁					
	A ₂					
	B					
PANOSCREEN LOT:	I					
	II					
	III					
DAT	Poly					
	IgG					
	C3					