

KPNW PPMP OBGYN Compe	tency Assessment 202	3 PERMANENTE
Employee Name (print):AliSon	Uyemura, ms	
Signature of Employee: Clin Ugenn		
Employee MUST also complete Kn	owledge Test: <u>www.medtraining.o</u> password=changeme	<u>rq</u> → login=email@kp.org,
	ach step of the procedure. Is with a question or observe the emp of as satisfactory. Ilowed, then mark unsatisfactory and p satisfactory rating. Describe actions to vation below on the 'Direct Observation'. In any one copy via endone copy	ployee resolving a problem.  perform necessary action during aken in comments section.  on/Assessment' line.
Microscopy Procedure (cross out if NOT Perform	ned): FERNING	T&Y
Assessment of Specimen Handling -proper specimen collection -specimen and SLIDE are tabeled appropriately -specimen properly prepared for testing (centrifuged/in saline etc.)	V∑ Salisfactory □Unsatisfactory	Satisfactory Unsatisfactory
Assessment of Test performance -microscope properly adjusted and clean -known result obtained	Satisfactory Unsatisfactory	Satisfactory Unsatisfactory
Assessment of Quality Control/Maintenance -verify proper storage of reagents and supplies -reagents in date/labeled appropriately -microscope maintenance complete	☑ Satisfactory □Unsatisfactory	Satisfactory Unsatisfactory
Assessment of Result entry into KPHC -Correct patient selected -Result Section filled out completely	Satisfactory Dunsatisfactory	Satisfactory Unsatisfactory
Assessment of Problem Solving skills Troubleshooting issues with Microscope	Satisfactory Dunsatisfactory	Satisfactory □Unsatisfactory
Direct Observation/Assessment by: \(\frac{1}{2}\)	shin Livingston, n	<u>が</u> Date: <u>4/21/23</u>
Comments:		