Kaiser Permanente Medical Center, Transfusion Service, Laboratory, 2425 Geary Blvd., San Francisco CA 94115.

## VISION Maintenance and QC Log

													a																		
Month: Year:												VISION #1 SN 50003189																			
<b>NOTE:</b> $\sqrt{1}$ for each task perfe	_	ed.	Mor	nthly	Ma	inter	nanc	ce sł	noul	d be	per	forn	ned	first,	, the	n W	/eek	ly th	en [												
Procedure	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Daily Maintenance																															
Empty Card Waste & Clean																															
Empty Liquid Waste													ĺ																		
Refill Saline & DI H2O																															
Daily Probe Maintenance (Remove partial cards from Incubator when cleaning probe)																															
Change Dil 2 & Dil 2+																															
Discard & Change RBC Rgts (every 5 days and on expiration day)																															
Check cards & dilution trays on- board inventory																															
Daily Rgt QC Performed by:																															
Daily Rgt QC Reviewed by:																															
Restart Analyzer (at end of Swing shift)																															
Weekly Maintenance (Thursday	)																														
Weekly Liquid System, Decomtamination, Pump Test																															
Performed by:																															
Weekly Maintenance (Thursday	and	on o	day c	of exp	oirati	on)																									
Change ALBA Q Check and new DAT POS & DAT NEG controls																															
Alba Q Check Lot#																															
Expiration date																															
Performed by:																															
Monthly Maintenance																															
Instrument Cleaning																															
Image System Cleaning																															
Data Backup																															
Performed by:																															
COMMENTS:																															

Supervisor Review:\_\_\_\_

Date:\_\_\_\_