

Assay	Negative	Equivocal (2 out of 3)	Re-test in Duplicate Zone (2 out of 3)	Hot Zone Positive
HAV IgG	<=0.99	none	none	>=1.00
HAV IgM	<=0.79	0.80-1.19	0.80-1.19	>=1.20
HBc Total	<=0.79	none	0.80-1.19	>=1.00
HBc IgM	<=0.79	0.80-1.19	0.80-1.19	>=1.20
AUSAB	<=7.99	8.00-11.99	8.00-11.99	>=12.00
HCV	<=0.79	0.80-0.99	0.80-0.99	>=1.00

Assay		Not Confirmed	Re-test with Dilution Zone	Confirmed
HBsAg		<=0.99	> =1.00	>=1.00 and CONF
Conf	Neat	1. <0.70 and NO % Neut 2. <10 and <50% Neut	>=10 and <50% Neut	>=0.70 and >=50% Neut
Conf	Dilu 1:500	1. 25 ul patient sample + 475 ul HBsAg Qual Confirmatory Manual Dil 2. 20 ul of Dilution #1 + 480 uL of HBsAg Qual Confirmatory Manual Dil	<0.70 and NO % Neut	>=0.70 and >=50% Neut
Conf	Dilu 1:20000	25 ul of 1:500 Dilution + 975 ul HBsAg Qual Confirmatory Manual Dil	1. <0.70 and NO % Neut 2. >= 0.70 and <50% Neut	>=0.70 and >=50% Neut