## SERVICE UPDATE



CAMPUS AFFECTED: PI	roctor		
	P	AFFECTS PATIE	NT CARE
TEST NAME: Erythrocyte Sedimentation Rate (ESR)			DATE OF CHANGE: 3/5/2018
REPLACING: NA	T		
SPECIMEN TYPE: 5ml Lavender EDTA Tube		avender EDTA Tube	
ORDERING CODES:	EPIC:	LAB322	NEW: NA
	SQ:	ESR	NEW: NA
OTHER USEFUL			
INFORMATION:			

In an effort to transition to the regionalization of laboratory services of in the Peoria region, all Erythrocyte Sedimentation Rate (ESR) testing will now be performed at Methodist campus laboratory.

Literature shows that C-Reactive Protein (CRP) is a more sensitive and accurate reflection of acute phase inflammation compared to ESR. CRP is recommended and still offered at the Proctor campus (EPIC Code: LAB149 SQ CODE: CRP)

Specimen stability requirements are: Ambient: 8 hours; refrigerated: 12 hours.

Testing is performed Sunday through Saturday.

Reference: https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4653962/

For further Information Refer to UnityPoint Methodist Laboratory User Manual

For Technical Information Contact:

For Scientific and Clinical Information Contact

Proctor Lab (309) 691-1060

Lori Racsa, DO, Clinical Pathologist 309-671-2181

The Do

Ron Fitzgerald, MLS (ASCP) Laboratory Asst. Manager 309-683-

Kon Johnson

6151

Adam Quinn, DO Medical Director 309-691-1060

Dana A Spears, Regional Director of Laboratories 309-672-4930