M-PO-TS-0102-04

PREADMIT PATIENTS - BLOOD BANK SPECIMENS (EXTENDED EXPIRATION)

☒ St. Joseph Medical Center Tacoma, WA
☒ St. Francis Hospital Federal Way, WA

St. Clare Hospital Lakewood, WASt. Anthony Hospital Gig Harbor, WA

☐ St. Elizabeth Hospital Enumclaw, WA☐ Highline Medical Center☑ PSC

PRINCIPLE

To ensure that compatible blood is available at the time of surgery, patients scheduled for surgery may have their pretransfusion specimens drawn up to fourteen days before the surgery date, provided certain regulatory requirements are met.

RELATED DOCUMENTS

R-W-TS0123	Blood Band – Typenex Barcoded
R-W-TS0124	Blood Band – Typenex R3 Attachment
R-W-TS0125	Blood Bands – When to Remove Them
M-W-TS-0100	Labeling Pre-Transfusion Specimens
M-F-TS1082	Preadmit Questions for Blood Bank Specime

M-F-TS1082 Preadmit Questions for Blood Bank Specimens Form

POLICY

- 1. Answers to all questions will be documented on the Preadmit Questions for Blood Bank Specimens Form.
- At the time of draw, the person who collects the blood specimen (phlebotomist) must <u>ask</u> the patient or the patient's caregiver accompanying the patient certain questions:
 - Has the patient been transfused during the last 90 days?
 - Has the patient been pregnant during the last 90 days?
 - What is the date of the surgery, procedure, or transfusion?
- 3. The patient or patient's caregiver must state that the patient:
 - Has not been transfused in the last three months
 - Females must state that they have not been pregnant in the last three months.
- 4. If the conditions in #3 have been met, the patient may be drawn for Type & Screen or Type & Crossmatch for up to 14 days before the date of the surgery, procedure, or transfusion.
- 5. If the conditions in #3 are NOT met (i.e., the patient HAS been transfused or pregnant within the last 90 days) the patient <u>cannot</u> be drawn early. The sample MUST be collected within 1-3 days prior to the surgery, procedure, or transfusion date. For instance, drawing a patient on Friday is acceptable for a Monday surgery 3 days later.
- 6. The specimen and the questionnaire form will be sent together to the SJMC Blood Bank
- 7. The Transfusion Service will retain this documentation for five years.

REFERENCES

AABB Standards for Blood Banks and Transfusion Services, current edition

G:\Lab\LAB\Document Control\Transfusion Service Active\Intranet - Loaded\Preadmit Patients - Blood Bank Specimens (Extended Expiration)-04.doc	Effective Date: 12/02/2014	Page 1 of 2	
Unauthorized use or copying of this document is prohibited by FHS.			

DOCUMENT APPROVAL **Purpose of Document / Reason for Change:** 1. Updated process due to Epic and loss of manual requisitions where phlebotomists could document patient history. 2. Added related documents section 3. Added references 4. Updated format 5. Changed to a multi-facility document rather than a regional one 6. Changed title to make it easier to search. Original title was "Pre-Transfusion Specimens on Pre-Admit Patients" No significant change to process in above revision. Per CAP, this revision does not require further Medical Director approval. □ Date: 12/1/2014 Karie Wilkinson, MD 12/1/14 Committee **Medical Director Approval Approval** specific document which is used at Date (Electronic Signature) only one facility