

IRICELL 3000 iChem Velocity Maintenance Log

Serial # V01535

Month/Year

DAILY	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Empty and Clean Waste Container																															
Clean Sample Transport Module																															
Clean Instrument Surfaces																															
Check test Strip & Wash Reagent																															
Run Daily QC																															
Replace desiccant																															
Initials																															



3 X WEEKLY (M, W, F) (as needed)	Week 1 Date	Week 2 Date	Week 3 Date	Week 4 Date	Week 5 Date
Clean Strip Provider Module (3x weekly M, W, F)					
Initials					

WEEKLY	Week1	Week 2	Week3	Week 4	Week 5
Clean Wash Station Bath					
Clean Strip Conveyor System					
Clean SGCC Module					
Initials					

As Needed	Date	Date	Date	Date	Date	Date	Date	Date	Date	Date
Replace Wash Solution and Filter										
Clean Sample Tube Detector										
Clean Barcode Reader Window										
Clean Optical Sensors on Sample Transport Module										
Initials										

QUARTERLY	QTR Date
Perform Reflectance Check	
Run Liquid CalCheck	
Initials	

CORRECTIVE ACTION	
Problem	Date/Time/Initials
Corrective Action	