



# **Entering/Exiting the BSL III with Minimal Personal Protective Equipment**

# BSL III PAPR with minimal PPE

- Entering the BSL III room 104 gowning room
  - Scan your ID badge to enter room 103
  - Sign into the 'BSL<sub>3</sub> Suite Activity Log' located outside of the door to room 104
    - Record date, time in, name and check off the box for 'not working with infectious agent'

BSL3 SUITE ACTIVITY LOG

Entry Date	Time In	Time Out	Name (PRINT) RO Authorized Personnel	Name (PRINT) Non-authorized Personnel [Informational Handout for Visitors, Maintenance, Service and Housekeeping Personnel on File]	CHOOSE ONLY ONE EXPLANATION	
					Work with infectious agent or material (PPE Required) (✓)	Not working with infectious agent or material [e.g., housekeeping / Maintenance / Service / Repair] (✓)
Date 9 / 9 /2013	Time In 10:24	Time Out	E. St Germain			<input checked="" type="checkbox"/>
Date / /2013	Time In	Time Out				

- Scan your ID badge to enter room 104

# BSL III PAPR with minimal PPE

- PPE Required to enter BSL III
  - Shoe Covers
  - One pair of gloves
- Entering the BSL III room 105
  - Check the ping pong ball air sensor over the door to room 105.
    - If you do not see the ball, you may enter the room.
    - If you see the ball, do not enter and contact your supervisor.
  - Enter your 5 digit code and press \* to enter.



# BSL III PAPR with minimal PPE

- Exiting the BSL III
  - Leave the BSL III and step into the yellow tape zone in room 104.
  - Remove one shoe cover and step outside of the yellow zone. Repeat for the second shoe cover. Place in the biohazard box.
  - Remove gloves and place in the biohazard box.
  - Wash your hands.
  - Leave room 104 and sign out of the Activity Log in room 103 by entering the time out.
  - Exit room 103 by scanning your badge.

# BSL III Entry/Exit with minimal PPE

Please take the quiz online: Passing grade is 80%

## Questions?

Please send an email to:

Eric St Germain, BT Coordinator

Beth Kolp, RO, LSO

Kevin Sohner, ARO

Tammy Bannerman, ARO

Larry King, PI