

**Blood Donor Center Daily QC Sheet**  
**In-House or Mobile Drive (Circle One)**

HemoCue 201+ OI 44-210	ID Number	Clean Holder (√)	Level 1 (g/dL) LOW	Level 2 (g/dL) NORM	Level 3 (g/dL) HIGH	Collection Area Temp (15-30°C)
			Lot:	Lot:	Lot:	
			Assigned Exp:	Assigned Exp:	Assigned Exp:	
			Manufacturer's Exp:	Manufacturer's Exp:	Manufacturer's Exp:	
			Range:	Range:	Range:	
						Tech:
						Tech:
						Tech: °C
						Tech:

Genesis QC/ Verification OI 44-218	Weight	Ranges	ECN	ECN	ECN	ECN	ECN	ECN	ECN	ECN	Tech	
	500g	499-501g										
	1 <sup>st</sup> Complete Unit	705-799g										

  

Genesis QC/ Verification OI 44-218	Weight	Ranges	ECN	ECN	ECN	ECN	ECN	ECN	ECN	ECN	Tech
	500g	499-501g									
	1 <sup>st</sup> Complete Unit	705-799g									

Suretemp Plus Thermometer QC OI 44-209	ECN	Self-Test (Sat/Unsat)	Clean (Use √)	Tech	Blood Pressure Monitor QC OI 44-211	ECN	Self-Test (Sat/Unsat)	Clean (Use √)	Tech	

  

Item	Lot Number	Manufacturer Exp Date	New Exp Date	Tech	OHAUS Electronic Scale QC OI 44-007	ECN	100g (+/-1)	200g (+/-1)	500g (+/- 1)	Tech
Hemocue Cuvettes										
Hemocue Cuvettes										

**Notes:**  
 - Notify Technical Supervisor of any discrepancies in quality control, lot numbers or range changes.  
 -Sat-Satisfactory/ -Unsat-Unsatisfactory

Comments: \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_  
 \_\_\_\_\_

