

HFC-A (HFC) 2009



Advancing Excellence

Last Updated:

July 7, 2009



KIT 21743561 6 01 67

CAP # 2471701 - 01 SEQ # 01

Products: HFC

OLYMPIC MEDICAL CENTER

STEVEN C BLACKHAM MT

TEL# 1-360-4177778

FAX# 360-417-7646

Page 1

Results must be received at the CAP no later than midnight, Central Time by the due date below:

Online: www.cap.org (preferred method)

Fax: 866-FAX-2CAP (866-329-2227)

Hemocytometer Fluid Count Survey Result Form

Enter Your Results Online at www.cap.org in Four Easy Steps!

Step 1: Opt In

Your laboratory must establish a laboratory web account to submit results online.

Step 2: Request Account

Laboratory staff who will enter results online must first establish a personal web account.

Step 3: Request Access

Once a personal web account is established, laboratory staff can request access to their laboratory's information.

Step 4: Enter Data!

For help, view the online demo and user guide on the CAP website.

© CAP 2009

Hemocytometer Fluid Count

Method	Exception Code		
⁰¹⁰ <input type="radio"/> 1368 Cell-VU Counting Chamber (Millenium Sciences) <input type="radio"/> 1071 Hemocytometer - bright field <input type="radio"/> 1070 Hemocytometer - phase <input type="radio"/> 1115 Hycor Kova System <input type="radio"/> 0010 Other, specify in final section	⁰²⁰ <input type="radio"/> 11 <input type="radio"/> 22 <input type="radio"/> 33		
Procedure	HFC-01	HFC-02	HFC-03
RBC/ μ L (RBC/ mm^3)	⁰³⁰ <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> .	⁰⁴⁰ <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> .	⁰⁵⁰ <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> .
WBC/ μ L (WBC/ mm^3)	⁰⁶⁰ <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> .	⁰⁷⁰ <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> .	⁰⁸⁰ <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> <input type="text"/> .

4426

